

Product Datasheet

SMAD2 (Phospho-S465/S467) antibody (orb639382)

Description	Rabbit polyclonal antibody to SMAD2 (Phospho-S465/S467)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S465/S467 of human Phospho-Smad2-S465/S467
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	60kDa
Uniprot ID	Q15796
Entrez	4087
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.