

Product Datasheet

RPA2 (Phospho-T21) antibody (orb639384)

Description	Rabbit polyclonal antibody to RPA2 (Phospho-T21)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding T21 of human RPA2.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P15927
Entrez	6118
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.