

Product Datasheet

PRKD1 (Phospho-S916) antibody (orb639415)

Description	Rabbit polyclonal antibody to PRKD1 (Phospho-S916)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S916 of human PKC mu(NP_002733.2?).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	102kDa
Uniprot ID	Q15139
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.