

Product Datasheet

Tau (Phospho-Ser416) antibody (orb639419)

Description	Rabbit polyclonal antibody to Tau (Phospho-Ser416)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A phospho specific peptide corresponding to residues surrounding Ser416 of human Tau.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	95kDa
Uniprot ID	P10636
Entrez	4137
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.