

Product Datasheet

RPS6KA1 (Phospho-T573) antibody (orb639421)

Description	Rabbit polyclonal antibody to RPS6KA1 (Phospho-T573)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T573 of human p90Rsk/RSK1/RPS6KA1(NP_002944.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	72kDa/81kDa/82kDa/83kDa
Uniprot ID	Q15418
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.