

Product Datasheet

RPS6KB1 (Phospho-T421/S424) antibody (orb639424)

Description	Rabbit polyclonal antibody to RPS6KB1 (Phospho-T421/S424)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T421 & S424 of human RPS6KB1(NP_001258989.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70kDa
Uniprot ID	P23443
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.