

Product Datasheet

IKBKG (Phospho-S376) antibody (orb639429)

Description	Rabbit polyclonal antibody to IKBKG (Phospho-S376)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S376 of human IKK γ (NP_003630.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	65kDa
Uniprot ID	Q9Y6K9
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.