

Product Datasheet

CD105 Antibody / Endoglin (orb639500)

Description

CD105/Endoglin is a Type I membrane

part of the TGF-beta...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A partial recombinant protein corresponding to amino acids 74-251 of human CD105/Endoglin was used as the immunogen for the Endoglin antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Endoglin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Titering of the Endoglin antibody may be required for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P17813

Hazard

Information

This Endoglin antibody is available for research use only.

Dilution Range

Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.