

Product Datasheet

PD-L1 antibody (PE) (orb639520)

Description	Mouse monoclonal antibody against PD-L1 conjugated to PE
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	PE
Tested Applications	FACS, IF
Immunogen	A portion of amino acids 39-191 from the human protein was used as the immunogen for this PD-L1 antibody.
Preservatives	500 ul at 0.1 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the PD-L1 antibody at 2-8°C, protected from light.
Note	For research use only
Application notes	Optimal dilution of the PD-L1 antibody should be determined by the researcher.
Formula	500 ul at 0.1 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	Q9NZQ7
Hazard Information	This PD-L1 antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells in 0.1ml, Immunofluorescence: 1-2ug/ml
Expiration Date	12 months from date of receipt.