

Product Datasheet

Clathrin Heavy Chain Antibody (orb639592)

Description

Recognizes protein of 192kDa, which is

is composed of t...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IF

Applications
Immunogen

Recombinant full-length human Clathrin Heavy Chain protein was used as the immunogen for this Clathrin Heavy Chain antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the Clathrin Heavy Chain antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

Q00610, **P53675**

Hazard
Information

This Clathrin Heavy Chain antibody is available for research use only.

Dilution Range

Immunofluorescence: 1-2ug/ml

Expiration Date

12 months from date of receipt.