


Product Datasheet

HCMV antibody (orb6396)

Description	HCMV antibody
Species/Host	Rabbit
Reactivity	HCMV
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	recombinant HCMV protein
Target	Cytomegalovirus Glycoprotein B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	74 kDa
Uniprot ID	F5HB81
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lymphocytes.