

Product Datasheet

**p120 Catenin Antibody (pTyr96) / CTNND1 /
Catenin delta 1 (orb639606)**

Description

The membrane associated protein

p120cas) was identified as a tyro...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FACS, IF

Applications
Immunogen

The amino acid sequence surrounding phosphorylated Tyr96 was used as the immunogen for the p120 Catenin antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the p120°Catenin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

O60716

Hazard
Information

This p120 Catenin antibody is available for research use only.

Dilution Range

Flow cytometry: 1-2ug/million cells in 0.1ml, Immunofluorescence: 1-2ug/ml

Expiration Date

12 months from date of receipt.