

Product Datasheet

BTLA Antibody / CD272 (orb639611)

Description	B and T lymphocyte attenuator (BTLA) is an immunoglobulin domain-containing glycoprotein whose expre...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant extracellular human CD272/BTLA protein was used as the immunogen for this BTLA antibody.
Preservatives	0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the BTLA antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the antibody should be determined by the researcher.
Formula	0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Hazard Information	This BTLA antibody is available for research use only.
Dilution Range	ELISA (order BSA-free format for coating)
Expiration Date	12 months from date of receipt.