

Product Datasheet

HER2 Antibody / ErbB2 (orb639623)

Description

This MAb is specific and shows minimal cross-

family. Recepto...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

A recombinant human HER2 protein fragment within amino acids 311-462 was used as the immunogen for the HER2 antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the HER2 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P04626

Hazard

Information

This HER2 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml

Expiration Date

12 months from date of receipt.