

Product Datasheet

ETS1 Antibody (orb639644)

Description

ETS1 proto-oncogene is an important

proliferation and...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications

Immunogen

A recombinant human partial protein (amino acids 137-230) was used as the immunogen for the ETS1 antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the ETS1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P14921

Hazard

Information

This ETS1 antibody is available for research use only.

Dilution Range

ELISA (order BSA-free format for coating)

Expiration Date

12 months from date of receipt.