

Product Datasheet

GnRH Receptor Antibody / LHCG Receptor (orb639686)

Description

Recognizes an epitope on the

releasing hormone (GnRH) receptor ...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

Recombinant full-length human GNRHR protein was used as the immunogen for this recombinant GNRHR / LHCGR antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the recombinant GNRHR / LHCGR antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Rabbit IgG, kappa

Clonality

Recombinant

Purity

Protein A affinity chromatography

Uniprot ID

P22888, **P30968**

Hazard
Information

This recombinant GNRHR / LHCGR antibody is available for research use only.

Dilution Range

ELISA (order BSA-free format for coating)

Expiration Date

12 months from date of receipt.