

Product Datasheet

HLA-DRB1 Antibody (orb639709)

Description

This MAb reacts with a 28kDa chain of HLA-

molecules. It does ...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FACS

Applications

Immunogen

Human REH cells (a non-T, non-B leukemia cell line) were used as the immunogen for the HLA-DRB1 antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the HLA-DRB1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P01911

Hazard

Information

This HLA-DRB1 antibody is available for research use only.

Dilution Range

Flow cytometry: 1-2ug/million cells

Expiration Date

12 months from date of receipt.