

Product Datasheet

Gamma Catenin Antibody / Plakoglobin (orb639757)

Description

It recognizes a protein of 80-87kDa,

catenins (-T-catenin, is locat...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FACS, IF, WB

Applications

Immunogen

Recombinant human full length protein was used as the immunogen for the recombinant Catenin gamma antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the recombinant Plakoglobin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Recombinant

Purity

Protein G affinity chromatography

Uniprot ID

P14923

Hazard

Information

This recombinant Plakoglobin antibody is available for research use only.

Dilution Range

ELISA (order BSA-free format for coating),Flow cytometry: 1-2ug/million cells in 0.1ml,Immunofluorescence: 1-2ug/ml,Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.