

Product Datasheet

Cytokeratin 14 Antibody / CK14 / KRT14 (orb639775)

Description

Cytokeratin 14 (CK14) belongs to the

molecular weight kerat...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FACS, IHC-P

Applications

Immunogen

A recombinant human partial protein (amino acids 351-472) was used as the immunogen for this Cytokeratin 14 antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the Cytokeratin 14 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Hazard

Information

This Cytokeratin 14 antibody is available for research use only.

Dilution Range

Flow cytometry: 1-2ug/million cells in 0.1ml, Immunohistochemistry (FFPE): 1-2ug/ml

Expiration Date

12 months from date of receipt.