

Product Datasheet

LSP1 Antibody / Lymphocyte Specific Protein 1 (orb639785)

Description

LSP1 is a hematopoietic-expressed

leukocyte-specific (includin...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

A recombinant human partial protein (amino acids 174-278) was used as the immunogen for this LSP1 antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the LSP1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P33241

Hazard

Information

This LSP1 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml

Expiration Date

12 months from date of receipt.