

Product Datasheet

Uroplakin 1A Antibody / UPK1A (orb639793)

Description

The protein encoded by this gene is a

superfamily, also known as the t...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC-P

Applications

Immunogen

Amino acids 114-173 were used as the immunogen for the Uroplakin 1A antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the Uroplakin 1A antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

000322

Hazard

Information

This Uroplakin 1A antibody is available for research use only.

Dilution Range

ELISA (order BSA-free format for coating), Immunohistochemistry (FFPE): 1-2ug/ml

Expiration Date

12 months from date of receipt.