

Product Datasheet

MMP3 Antibody / Matrix Metalloproteinase 3 (orb639814)

Description

The matrix metalloproteinases (MMP)

responsible for the degradatio...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

Amino acids 317-327 from the human protein were used as the immunogen for this Matrix Metalloproteinase 3 antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the Matrix Metalloproteinase 3 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, lambda

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID
P08254
Hazard
Information

This Matrix Metalloproteinase 3 antibody is available for research use only.

Dilution Range

ELISA: 2-4ug/ml for coating (order BSA/azide-free format)

Expiration Date

12 months from date of receipt.