

Product Datasheet

c-Myc Antibody (orb639824)

Description

It recognizes a transcription factor of 64-67kDa,

aa ...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FACS

Applications
Immunogen

Amino acids 408-439 were used as the immunogen for the recombinant c-Myc antibody. The specific epitope has been localized to amino acids 410-419 (EQKLISEEDL).

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the recombinant c-Myc antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Rabbit IgG, kappa

Clonality

Recombinant

Purity

Protein A affinity chromatography

Uniprot ID

P01106

Hazard
Information

This recombinant c-Myc antibody is available for research use only.

Dilution Range

ELISA (order BSA-free format for coating), Flow cytometry: 1-2ug/million cells in 0.1ml

Expiration Date

12 months from date of receipt.