

Product Datasheet

NME1 Antibody / NM23 (orb639833)

Description

The nm23 gene, a potential suppressor of

differential hybri...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A recombinant full-length human NME1 protein was used as the immunogen for this NME1 antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the NME1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P15531

Hazard

Information

This NME1 antibody is available for research use only.

Dilution Range

Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.