

Product Datasheet

Cytochrome C Antibody (orb639865)

Description

In mammalian cells, the highly conserved

the mitochond...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FACS, IF, IHC-P, WB

Applications
Immunogen

A recombinant full-length human protein was used as the immunogen for this recombinant Cytochrome C antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the recombinant Cytochrome C antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Rabbit IgG

Clonality

Recombinant

Purity

Protein A affinity chromatography

Uniprot ID

P99999

Hazard
Information

This recombinant Cytochrome C antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 0.25-0.5ug/ml, Flow cytometry: 1-2ug/million cells, Western blot: 1-2ug/ml, Immunofluorescence: 1-2ug/ml

Expiration Date

12 months from date of receipt.