

## Product Datasheet

### Bcl-X Antibody (orb639905)

**Description**

Recognizes a protein of 27kDa, identified as the

wit...

**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

ELISA

**Applications**
**Immunogen**

A recombinant human partial protein (amino acids 1-135) was used as the immunogen for this Bcl-X antibody.

**Preservatives**

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

**Storage**

Store the Bcl-X antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

**Note**

For research use only

**Application notes**

Optimal dilution of the antibody should be determined by the researcher.

**Formula**

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

**Isotype**

Mouse IgG2b, kappa

**Clonality**

Monoclonal

**Purity**

Protein G affinity chromatography

**Uniprot ID**

**Q07817**

**Hazard**
**Information**

This Bcl-X antibody is available for research use only.

**Dilution Range**

ELISA (order BSA-free format for coating)

**Expiration Date**

12 months from date of receipt.