

Product Datasheet

L-Selectin Antibody / CD62L (orb639920)

Description

Selectins, also designated CD62 antigens,

proteins involve...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

Supernatant from phorbol myristic acid activated human peripheral blood leukocytes was used as the immunogen for the L-Selectin antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the CD62L antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P14151

Hazard
Information

This CD62L antibody is available for research use only.

Dilution Range

ELISA (order BSA-free format for coating)

Expiration Date

12 months from date of receipt.