

Product Datasheet

BRCA1 Antibody (Breast Marker) (orb639942)

Description

This gene encodes a nuclear

maintaining genomic stability, and i...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications

Immunogen

A recombinant human partial protein (amino acids 445-620) was used as the immunogen for this BRCA1 antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the BRCA1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P38398

Hazard

Information

This BRCA1 antibody is available for research use only.

Dilution Range

ELISA (coating): 2-4ug/ml (order BSA-free format)

Expiration Date

12 months from date of receipt.