

Product Datasheet

VLDLR Antibody / VLDL Receptor (orb639995)

Description

VLDLR (very low density lipoprotein

receptor gene family, which inc...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

A recombinant protein fragment from the C-terminal region of VLDLR was used as the immunogen for this recombinant VLDL Receptor antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the recombinant VLDL Receptor antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Rabbit IgG, kappa

Clonality

Recombinant

Purity

Protein A affinity chromatography

Uniprot ID
P98155
Hazard
Information

This recombinant VLDL Receptor antibody is available for research use only.

Dilution Range

ELISA (order BSA-free format for coating)

Expiration Date

12 months from date of receipt.