

Product Datasheet

GITR Antibody / TNFRSF18 (orb640014)

Description

The tumor necrosis factor receptor (TNFR)

type I transmem...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

Human full-length TNFRSF18 protein was used as the immunogen for this TNFRSF18 antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the TNFRSF18 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG3, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID
Q9Y5U5
Hazard
Information

This TNFRSF18 antibody is available for research use only.

Dilution Range

ELISA (order BSA-free format for coating & Azide)

Expiration Date

12 months from date of receipt.