

Product Datasheet

CD80 Antibody (orb640045)

Description

T cell proliferation and lymphokine production

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

A recombinant human partial protein (amino acids 35-142) was used as the immunogen for this CD80 antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the CD80 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID
P33681
Hazard
Information

This CD80 antibody is available for research use only.

Dilution Range

ELISA (order BSA-free format for coating)

Expiration Date

12 months from date of receipt.