

Product Datasheet

Clostridium difficile Toxin A Antibody (orb640085)

Description

This mAb reacts with C. difficile Toxin A,

toxin, Pseudomona...

Species/Host

Mouse

Reactivity

Bacteria

Conjugation

Unconjugated

Tested

ELISA, IF

Applications
Immunogen

C. difficile Toxin A was used as the immunogen for the Clostridium difficile Toxin A antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the Clostridium difficile Toxin A antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P16154

Hazard
Information

This Clostridium difficile Toxin A antibody is available for research use only.

Dilution Range

ELISA (order BSA-free format for coating), Immunofluorescence: 1-2ug/ml

Expiration Date

12 months from date of receipt.