

## Product Datasheet

### c-RET Antibody (orb640098)

## Description

The Ret proto-oncogene is structurally related to

transmembran...

## Species/Host

Mouse

## Reactivity

Human

## Conjugation

Unconjugated

## Tested

ELISA

## Applications

## Immunogen

A recombinant human partial protein (amino acids 702-848) was used as the immunogen for this c-RET antibody.

## Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

## Storage

Store the c-RET antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

## Note

For research use only

## Application notes

Optimal dilution of the antibody should be determined by the researcher.

## Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

## Isotype

Mouse IgG1, kappa

## Clonality

Monoclonal

## Purity

Protein G affinity chromatography

## Uniprot ID

**P07949**

## Hazard

## Information

This c-RET antibody is available for research use only.

## Dilution Range

ELISA (order BSA-free format for coating)

## Expiration Date

12 months from date of receipt.