

Product Datasheet

Interferon gamma Antibody / IFNG (orb640107)

Description

Recognizes a protein of 20-25kDa,

is specific to human IFN-...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

Purified recombinant human IFNG was used as the immunogen for the recombinant Interferon gamma antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the recombinant Interferon gamma antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Rabbit IgG, kappa

Clonality

Recombinant

Purity

Protein A affinity chromatography

Uniprot ID

P01579

Hazard
Information

This recombinant Interferon gamma antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml

Expiration Date

12 months from date of receipt.