

Product Datasheet

OX40 Antibody / CD134 (orb640139)

Description

CD134 is a type I integral membrane protein.

Species/Host

Rat

Reactivity

Mouse

Conjugation

Unconjugated

Tested Applications

FACS, IF

Immunogen

Purified recombinant mouse OX40 antigen was used as the immunogen for the CD134 antibody.

Preservatives

0.2 mg/ml with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the CD134 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the antibody should be determined by the researcher.

Formula

0.2 mg/ml with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Rat IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P47741

Hazard Information

This CD134 antibody is available for research use only.

Dilution Range

This antibody does not block binding of OX40L, Flow cytometry: 1-2ug/million cells in 0.1ml, Immunofluorescence: 1-2ug/ml

Expiration Date

12 months from date of receipt.