

Product Datasheet

Human Novel Coronavirus Nucleoprotein protein (orb640276)

Description

Recombinant human Novel Coronavirus Nucleoprotein(N) protein

Endotoxins

Not test.

Preservatives

Lyophilized from a 0.2 µm filtered 20 mM Tris-HCl, 0.5 M NaCl, 6% Trehalose, pH 8.0

Form/Appearance

Lyophilized powder

Storage

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Tag

N-terminal 6xHis-sumostar-tagged

Note

For research use only

Application notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Protein Sequence

MSDNGPQNQRNAPRITFGGPSDSTGSNQNGERSGARSQRRPQGLPNNTASW
FTALTQHGKEDLKFRGQGVPIINTNSSPDDQIGYYRRATRRIRGGDGKMKDLSP
RWYFYLLGTGPEAGLPYGANKDGIWVATEGALNTPKDHIGTRNPANNAIVLQL
PQGTTLPKGFYAEGSRGGSQASSRSSRSRNSSRNSTPGSSRGTSPPARMAGNG
GDAALALLLLDRLNQLESKMSGKGQQQGGQTVTKKSAEASKKPRQKRTATKA
YNVTQAFGRRGPEQTQGNFGDQELIRQGTQDYKHWPQIAQFAPSASAFFGMSRI
GMEVTPSGTWLTYAAIKLDDKDPNFKDQVILLNKHIDAYKTFPPTEPKDKKKK
ADETQALPQRQKKQTVTLLPAADLDDFSKQLQSMSSADSTQA

Purity

Greater than 90% as determined by SDS-PAGE.

Source

Yeast

MW

67.6 kDa

Uniprot ID
P0DTC9
Expiration Date

6 months from date of receipt.