

Product Datasheet

PI3K p55 Gamma (Y463) antibody (orb640786)

Description	Rabbit polyclonal antibody to PI3K p55 Gamma (Y463)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide, corresponding to amino acids 430-480 of Human PI3K p85 α .
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 55 kDa
Uniprot ID	P27986, Q92569
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.