

Product Datasheet

DDDDK-tag (3B9) antibody (orb640801)

Description	Mouse monoclonal antibody to DDDDK-tag (3B9)
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide DDDDK-tag.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -22°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: 1:5000-1:10000
Expiration Date	12 months from date of receipt.