

Product Datasheet

Smad2/3 (Phospho-T8) antibody (orb640807)

Description	Rabbit polyclonal antibody to Smad2/3 (Phospho-T8)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human Smad2/3 around the phosphorylation site of Threonine 8.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 52 kDa
Uniprot ID	P84022, Q15796
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.