

Product Datasheet

DDDDK-tag antibody (orb640817)

Description	Rabbit polyclonal antibody to DDDDK-tag
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Synthetic peptide DDDDK-tag
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:5000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.