

Product Datasheet

DDDDK-tag (1A8) antibody (orb640818)

Description	Mouse monoclonal antibody to DDDDK-tag (1A8)
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Synthetic peptide DDDDK-tag.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: 1:5000-1:50000IP: 1:100-1:1000, IF: 1:100-1:1000
Expiration Date	12 months from date of receipt.