

Product Datasheet

DDDDK-tag antibody (orb640820)

Description	Mouse monoclonal antibody to DDDDK-tag
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide DDDDK-tag.
Preservatives	1 mg/ml in Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Uniprot ID	N/A
Dilution Range	WB: 1:2000-1:5000
Expiration Date	12 months from date of receipt.