

Product Datasheet

PARP antibody (orb640874)

Description	Rabbit polyclonal antibody to PARP
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant full length Human PARP.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 100, 110 kDa
Uniprot ID	P09874
Dilution Range	WB1:1000 - 1:3000, IHC-P: -P1:50 - 1:200, IF1:50 - 1:200IP1:50 - 1:200, ChIP1:50 - 1:200
Expiration Date	12 months from date of receipt.