

Product Datasheet

Aminoallyl-UTP (Cy3) (orb64091)

Description	Aminoallyl-UTP conjugated to Cy3.
Form/Appearance	solution in water; Color: pink; pH: 7.5 ± 0.5
Concentration	1.0 mM-1.1 mM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 550 nm, λ_{em} 570 nm, ϵ 150.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
Formula	C ₄₃ H ₅₆ N ₅ O ₂₂ P ₃ S ₂
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 1151.97 g/mol (free acid); Detected MW: 1151.21 g/mol (free acid)
SMILES	<chem>P(=O)(OP(=O)(O)O)(O)OP(=O)(O)OC[C@H]1O[C@H]([C@@H]([C@@H]1O)O)n1c(=O)[nH]c(=O)c(c1)/C=C/CNC(=O)CCCCCN1c2ccc(S(=O)(=O)O)cc2C(C)(C)/C1=C\C=C\C1=[N+](c2ccc(cc2C1(C)C)S(=O)(=O)[O-])CC</chem>
Expiration Date	12 months from date of receipt.