

Product Datasheet

CLUH antibody (orb640971)

Description	Rabbit polyclonal antibody to CLUH
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant protein of human CLUH.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 180 kDa
Uniprot ID	O75153
Dilution Range	WB1:1000 - 1:4000, IHC-P: -P1:50 - 1:200, IF1:50 - 1:200
Expiration Date	12 months from date of receipt.