

Product Datasheet

ZC3H14 antibody (orb641028)

Description	Rabbit polyclonal antibody to ZC3H14
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant protein of human ZC3H14
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 105 kDa
Uniprot ID	Q6PJT7
Dilution Range	WB 1:500 - 1:2000 , IHC-P: -P 1:50 - 1:200
Expiration Date	12 months from date of receipt.