

Product Datasheet

P4HA1 antibody (orb641085)

Description	Rabbit polyclonal antibody to P4HA1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 235-534 of human P4HA1 (NP_001017962.1).
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 61 kDa
Uniprot ID	P13674
Dilution Range	WB 1:1000 - 1:2000
Expiration Date	12 months from date of receipt.