

Product Datasheet

CTCF antibody (orb641323)

Description	Rabbit polyclonal antibody to CTCF
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, IP, WB
Immunogen	Recombinant full length Human CTCF.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 140 kDa
Uniprot ID	P49711
Dilution Range	WB 1:1000 - 1:2000, IHC-P: -P/, IF 1:50 - 1:200IP/, ChIP 1:50 - 1:100
Expiration Date	12 months from date of receipt.