

Product Datasheet

Aminoallyl-dUTP (ATTO-Rho101) (orb64145)

Description	Aminoallyl-dUTP conjugated to ATTO-Rho101.
Form/Appearance	solution in water; pH: 7.5 ± 0.5
Concentration	1.0 mM-1.1 mM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 587 nm, λ_{em} 609 nm, ϵ 120.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
Formula	C ₄₉ H ₅₇ N ₆ O ₁₇ P ₃
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 1094.94 g/mol (free acid); Detected MW: 1094.30 g/mol (free acid)
SMILES	<chem>P(=O)(OP(=O)(O)O)(O)OP(=O)(O)OC[C@H]1O[C@H](C[C@@H]1O)n1c(=O)[nH]c(=O)c(c1)/C=C/CNC(=O)CCCN(C(=O)c1c(C2=c3cc4c5=[N+](C(CC4)CCCc5c3O)c3c2cc2c4c3CCCN4CCC2)cccc1)C</chem>
Expiration Date	12 months from date of receipt.