

Product Datasheet

Aminoallyl-dUTP (ATTO-665) (orb64151)

Description	Aminoallyl-dUTP conjugated to ATTO-665.
Form/Appearance	solution in water; Color: blue; pH: 7.5 ± 0.5
Concentration	0.50 mM-0.55 mM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 662 nm, λ_{em} 680 nm, ϵ 160.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
Formula	C ₁₂ H ₂₀ N ₃ O ₁₄ P ₃ -ATTO 665
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 1127.22 g/mol (free acid)
SMILES	<chem>P(=O)(OP(=O)(O)O)(O)OP(=O)(O)OC[C@H]1O[C@H](C[C@@H]1O)n1c(=O)[nH]c(=O)c(c1)/C=C/CN*</chem>
Expiration Date	12 months from date of receipt.